## ClinCal<sup>®</sup> - Calibrator

**Plasma Calibrator lyophilised, for Antiasthmatics /**

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The mean values have been assigned with the ClinRep<sup>®</sup> Complete Kit for Antiasthmatics.

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### Intended use:

ClinCal<sup>®</sup> Plasma Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

### Reconstitution:

Add exactly 3.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

### Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma calibrator is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 24 hours when stored at 15 - 30 °C
  - at least 7 days when stored at 2 - 8 °C
  - at least 3 months when stored below -18 °C
- (avoid repeated freezing and thawing)

The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.98 %.

### Pack size:

ClinCal<sup>®</sup> Plasma Calibrator for Antiasthmatics: 1 x 3 ml, **order no.: 12013**

### Caution:

The human plasma which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the plasma calibrator should be considered as potentially infectious and treated with appropriate care.

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