IDK[®] Casomorphin/Gliadorphin- Peptide LC-MS/MS

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
IDK[®] Casomorphin/Gliadorphin peptides LC-MS/MS

For the determination in stabilised urine

Valid from 2016-12-01

KM8000

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	19
2. INTRODUCTION	19
3. MATERIAL SUPPLIED	20
4. MATERIAL REQUIRED BUT NOT SUPPLIED	20
5. STORAGE AND PREPARATION OF REAGENTS	21
6. STORAGE AND PREPARATION OF SAMPLES	21
7. ASSAY PROCEDURE	21
<i>Preparation of the cartridges</i>	22
<i>Test procedure</i>	22
8. CHROMATOGRAPHIC CONDITIONS: UPLC METHOD	22
9. MS/MS METHOD	23
<i>MRM transitions (m/z)</i>	23
<i>Calculation</i>	24
10. CHROMATOGRAM EXAMPLES	25
11. QUALITY CONTROL	28
<i>Reference range</i>	28
12. PERFORMANCE CHARACTERISTICS	28
<i>Precision and reproducibility</i>	28
<i>Sensitivity</i>	30
13. PRECAUTIONS	30
14. TECHNICAL HINTS	30
15. DISPOSAL	30
16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	31

1. INTENDED USE

The described ImmuTube® LC-MS/MS application is intended for the quantitative determination of casomorphin/gliadorphin peptides in stabilised urine. For *in vitro* diagnostic use only.

2. INTRODUCTION

Exorphins are bioactive peptides which can react with opiate receptors. These peptides are derived from e.g. the incomplete hydrolysis of the milk protein casein and the wheat protein gluten and gliadin. Casomorphin (BCM7 from casein) and gliadorphin-7 (from gluten) have been widely studied and are linked to gastrointestinal, neurological, and neuro-developmental disorders.

These short peptides are usually broken down by the enzyme dipeptidyl-peptidase IV (DPPIV). In people with autism, schizophrenia and other neurological disorders, either a genetic deficiency or inhibition of this enzyme, e.g. by drugs treating diabetes type 2, mercury or yeast, can occur.

Both ways, this proteolytic enzyme can't act properly and leaves incompletely digested peptides behind which are able to pass the blood-brain barrier after the uptake from the intestinal lining into the bloodstream. In the brain, these bioactive peptides can bind to opiate receptors and mimic the effects of opiate drugs like heroine and morphine. This can result in change of behavior, lack of focus and attention, sleepiness and even aggression and self-abuse.

Wheat and milk allergies are associated casomorphin/gliadorphin problems. Even if no milk or wheat allergies occur some people can still react negatively to these peptides.

Children with autism and ADHS can benefit from an elimination diet – gluten-free, casein-free (GfCf) diet. Physical and mental health can be improved.

This test for casomorphin/gliadorphin peptides is superior to any other tests because no nutritional regimen, e.g. abandonment of soy products, prior to sample collection is required.

3. MATERIAL SUPPLIED

Cat. No.	Label	Kit components	Quantity
KM8000LMA KM8000LMB	MOPHA A MOPHA B	Mobile Phase	1000 ml 1000 ml
KM8000KA	CAL 1 CAL 2	Calibrator 1 and 2, (2.2 ml lyophilised, for concentration see product specification)	5 vials each
KM8000KO	CTRL 1 CTRL 2	Control 1 and 2 (2.2 ml lyophilised, for concentration see product specification)	5 vials each
KM8000IS	INT STD	Internal Standard	1.5 ml
KM8000AC	ACTSOL	Activation reagent	2.5 ml
KM8000WL	WASHSOL	Wash solution	450 ml
KM8000SO	SOL A	Solution A	22 ml
KM8000EL	ELUSOL	Elution solution	55 ml
KM8000RE	RECSOL	Reconstitution solution	4 x 15 ml

For the preparation of the casomorphin/gliadorphin peptides LC-MS/MS application, tuning is necessary to find the ideal LC-MS/MS settings. The tuning solution (KM8000TU [5 ml], KM8000TS [2 ml]), the UPLC column (KM8000RP), the SPE cartridges (KM8000CK) and all other individual components can be ordered separately. Please ask for our single component price list.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- SPE cartridges 30 mg/1ml (Strata-X 33u from Phenomenex; KM8000CK)
- Standard laboratory glass vials, suitable for LC-MS/MS
- Calibrated precision pipettors and 10–1000 µl tips
- Solid phase extraction unit
- Evaporation unit (e.g. vacuum centrifuge or nitrogen distributor)
- Vortex
- LC-MS/MS equipment
- RP-C18 column, e.g. Aeris PEPTIDE 3,6u XB (15,0 × 2,1 mm) from Phenomenex (KM8000RP)
- Methanol p.a.
- Ultra pure water*

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5. STORAGE AND PREPARATION OF REAGENTS

- Lyophilised calibrators (**CAL 1** and **CAL 2**) and lyophilised controls (**CTRL 1** and **CTRL 2**) are stable at **-20 °C** until the expiry date stated on the label. Before use, the lyophilised calibrators and controls needed have to be reconstituted with 2.2 ml reconstitution solution (RECSOL). The concentration of casomorphin/gliadorphin-peptides slightly changes from lot to lot. The exact concentration is stated on the specification data sheet. **Reconstituted calibrators and controls are not stable and cannot be stored.**
- The **mobile phases** (MOPHA A, MOPHA B) and **solution A (SOL A)** have to be spiked with **0.1 % activation reagent (ACTSOL)** before use, e. g.

500 ml MOPHA + 500 µl ACTSOL

11 ml SOL A + 11 µl ACTSOL

Mobile phases and solution A spiked with activation reagent can be used for **4 weeks**. Prepare only the appropriate amount necessary for each run.

Attention: The **activation reagent** must be added under the **fume hood**. All vials to be used must be absolutely clean, detergent-free and preferably made of glass suitable for LC-MS/MS.

- All other test reagents are stable until the expiry date (see label of test package) when stored at **2–8 °C**.

6. STORAGE AND PREPARATION OF SAMPLES

In unstabilised urine, the concentration of casomorphin/gliadorphin peptides declines quickly, therefore only **stabilised urine** is suited as sample. Immundiagnostik offers urine sample stabilising tubes (order no DZ9030UT) specifically designed for this purpose. Samples stabilised using those tubes are stable for 48 h at room temperature. For longer storage, the sample must be frozen at < -17 °C.

7. ASSAY PROCEDURE

Control samples should be analysed with each run.

Bring all reagents and samples to room temperature (18–24 °C) and mix well.

Preparation of the cartridges

Conditioning of the cartridges with 1 ml methanol, followed by equilibration with 1 ml ultra pure water.

Test procedure

1.	Add to 1 ml SAMPLE, CAL 1, CAL 2, CTRL 1, CTRL 2 (stabilised urine sample, calibrators and controls) each 10 µl INT STD (internal standard), mix, pipette onto the cartridges and evacuate.
2.	Suck 4 x 1 ml wash solution (WASHSOL) through the column.
3.	Eluate with 500 µl elution solution (ELUSOL) into glass tubes.
4.	Dry the eluate , e. g. in a vacuum centrifuge.
5.	Reconstitute the dried sample with 200 µl activated SOL A .
6.	Inject 15 µl into the LC-MS/MS system.

8. CHROMATOGRAPHIC CONDITIONS: UPLC METHOD

Column material: e. g. Aeris PEPTIDE 3,6 u XB from Phenomenex

Column dimension: 2.1 × 150 mm

Flow rate: 0.2 ml/min

Column temperature: 30°C

Injection volume: 15 µl

Running time: 15 min

Gradient:

Time (min)	% A	% B
0	95	5
2	95	5
10	50	50
10,1	5	95
12	5	95
12,1	95	5
15	95	5

It is recommended to use a guard column/filter to extend the column's life. The guard column cartridge (KM8000VS) and guard column holder (KM8000VH) can be ordered separately from Immundiagnostik.

After the analysis, the separation column should be washed with ~ 20 ml 50% methanol. The column can be stored in 50% methanol.

9. MS/MS METHOD

(Exemplified for a Waters TQ MS tandem mass spectrometer)

Mode:	MRM
Polarity:	ESI ⁺
Capillary voltage (kV):	3
Cone voltage (V):	var.
Source temperature (°C):	150
Desolvation temperature (°C):	600
Gas flow at cone (l/h):	100
Desolvation gas flow (l/h):	800
Collision gas flow (ml/min):	0.15

We recommend to set the preferences for the scantime to result in 8–15 data points per peak.

MRM transitions (m/z)

β-Casomorphin 7 (molecular weight: 789.95 g/mol)

790.43 > 70.05	cone voltage: 52	collision energy: 70
790.43 > 136.1	cone voltage: 52	collision energy: 70

INT STD (β-Casomorphin 7 heavy labelled, molecular weight: 795.95 g/mol)

796.63 > 70.36	cone voltage: 45	collision energy: 60
796.63 > 74.95	cone voltage: 45	collision energy: 60

Casomorphin 1–3 (molecular weight: 425.50 g/mol)

426.27 > 70.09	cone voltage: 29	collision energy: 29
426.27 > 136.14	cone voltage: 29	collision energy: 31
426.27 > 263.22	cone voltage: 29	collision energy: 19

Casomorphin 1–4 (molecular weight: 522.62 g/mol)

523.13 > 70.08	cone voltage: 35	collision energy: 40
523.13 > 116.07	cone voltage: 35	collision energy: 40
523.13 > 136.10	cone voltage: 35	collision energy: 40

Casomorphin 1–4amid (molecular weight: 521.63 g/mol)

522.27 > 70.05	cone voltage: 32	collision energy: 39
522.27 > 136.06	cone voltage: 32	collision energy: 39

Exorphin A5 (molecular weight: 599.65 g/mol)

600.41 > 70.09	cone voltage: 23	collision energy: 43
600.41 > 136.13	cone voltage: 23	collision energy: 29
600.41 > 217.17	cone voltage: 23	collision energy: 17

Exorphin B5 (molecular weight: 594.67 g/mol)

595.35 > 132.04	cone voltage: 31	collision energy: 55
595.35 > 136.14	cone voltage: 31	collision energy: 50

Exorphin C (molecular weight: 591.72 g/mol)

592.32 > 70.05	cone voltage: 40	collision energy: 45
592.32 > 261.12	cone voltage: 40	collision energy: 40

Gliadorphin (molecular weight: 876.00 g/mol)

876.52 > 70.05	cone voltage: 55	collision energy: 61
876.52 > 115.17	cone voltage: 55	collision energy: 65
876.52 > 136.10	cone voltage: 55	collision energy: 60

Calculation


Linear regression can be used as model for evaluation of the results (in consideration of the internal standard). The two calibrator concentration points are connected by a straight line. The samples can be calculated using this obtained line.

10. CHROMATOGRAM EXAMPLES


β -Casomorphin 7 (molecular weight: 789.95 g/mol)

Sample Name: Casoglia150616_001

F7:MRM of 2 channels,ES+
790.43 > 136.1
5.926e+004


F7:MRM of 2 channels,ES+
790.43 > 70.05
2.175e+005


INT STD (β -Casomorphin 7 heavy labelled, molecular weight: 795.95 g/mol)

Sample Name: Casoglia150616_001

F8:MRM of 2 channels,ES+
796.63 > 74.95
1.517e+005


F8:MRM of 2 channels,ES+
796.63 > 70.36
2.215e+004


Casomorphin 1-3 (molecular weight: 425.50 g/mol)

Sample Name: Casoglia150616_001

F1:MRM of 2 channels,ES+
426.82 > 136.1
9.308e+003


F1:MRM of 2 channels,ES+
426.82 > 70.32
6.050e+003


Casomorphin 1-4 (molecular weight: 522.62 g/mol)**Sample Name: Casoglia150616_001**


F3:MRM of 3 channels,ES+
523.13 > 70.08
1.668e+005


F3:MRM of 3 channels,ES+
523.13 > 116.07
4.734e+004

**Casomorphin 1-4amid** (molecular weight: 521.63 g/mol)**Sample Name: Casoglia150616_001**


F2:MRM of 2 channels,ES+
522.27 > 136.07
1.409e+005


F2:MRM of 2 channels,ES+
522.27 > 70.05
2.280e+005







**Exorphin A5** (molecular weight: 599.65 g/mol)**Sample Name: Casoglia150616_001**

F6:MRM of 2 channels,ES+
600.06 > 70.05
6.321e+003

F6:MRM of 2 channels,ES+
600.06 > 136.07
1.126e+004

Exorphin B5 (molecular weight: 594.67 g/mol)**Sample Name: Casoglia150616_001****Exorphin C** (molecular weight: 591.72 g/mol)**Sample Name: Casoglia150616_001****Gliadorphin** (molecular weight: 876.00 g/mol)**Sample Name: Casoglia150616_001**

11. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Reference range

Based on Immundiagnostik studies of samples of apparently healthy persons (n = 18), the following mean values were estimated

Analyte	ng/ μ mol Creatinine
Casomorphin 1-3	1.406
Exorphin A5	2.033
β -Casomorphin 7	0.138
Gliadorphin	0.187
Casomorphin 1-4 amid	0.234
Casomorphin 1-4	0.188
Exorphin C	0.076
Exorphin B5	0.611

We recommend each laboratory to establish its own reference range.

12. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Inter-Assay (n = 10)

Analyte	Sample	ng/ml	CV [%]
Casomorphin 1-3	1	1.05	10.79
	2	4.55	2.41
Exorphin A	1	3.75	24.35
	2	11.37	9.27

Analyte	Sample	ng/ml	CV [%]
β-Casomorphin 7	1	1.56	12.32
	2	6.93	4.66
Gliadorphin	1	1.96	7.38
	2	7.80	7.14
Casomorphin 1-4 amid	1	1.44	11.40
	2	5.57	6.07
Casomorphin 1-4	1	2.09	5.13
	2	8.71	5.54
Exorphin C	1	2.84	5.16
	2	11.76	8.53
Exorphin B	1	2.88	5.37
	2	10.91	4.31

Intra-Assay (n = 4)

Analyte	Sample	ng/ml	CV [%]
Casomorphin 1-3	1	3.44	12.92
	2	1.14	3.63
Exorphin A	1	6.29	26.98
	2	3.05	39.92
β-Casomorphin 7	1	3.19	8.42
	2	1.52	11.88
Gliadorphin	1	3.17	8.03
	2	2.14	7.59
Casomorphin 1-4 amid	1	2.33	17.08
	2	1.34	15.58
Casomorphin 1-4	1	2.83	7.79
	2	2.06	7.77
Exorphin C	1	3.91	4.11
	2	3.05	4.69
Exorphin B	1	5.24	13.38
	2	2.72	8.84

Sensitivity

Detection limit

Analyte	pg/ml Urine
Casomorphin 1-3	1132
Exorphin A	1872
β -Casomorphin 7	420
Gliadorphin	520
Casomorphin 1-4 amid	624
Casomorphin 1-4	307
Exorphin C	473
Exorphin B	976

It should be noted that the determination of the detection limit depends not only on the application method but also on the instrument.

13. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Control samples should be analysed with each run.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

14. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date stated on kit label.
- The assay should always be performed according to the enclosed manual.

15. DISPOSAL

Mobile phase B, elution solution, solution A and activation reagent have to be disposed as non-halogenated solvents.

16. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK*[®] is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		