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# **IDK<sup>®</sup> IDO ELISA**

***Zur in-vitro-Bestimmung von Indolamin-2,3-Dioxygenase 1  
(IDO1)***

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### 14. LITERATUR

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; **8** (3): 289-95.

Cavia-Saiz M, Muñoz Rodríguez P, Llorente Ayala B, García-González M, Coma-Del Corral MJ, García Girón C. The role of plasma IDO activity as a diagnostic marker of patients with colorectal cancer. *Mol Biol Rep.* 2014 Apr; **41**(4):2275-9.

Ciorba MA: Indoleamine 2,3 dioxygenase in intestinal disease. *Curr Opin Gastroenterol.* 2013 Mar;**29**(2):146-52

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology.* 2013 Mar 1; **2**(3):e23428

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol.* 2012 March; **19**(3): 436–442

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**Manual**

# **IDK® IDO ELISA**


***For the in vitro determination of indoleamine 2,3-dioxygenase 1  
(IDO1)***

***For research use only***

Valid from 2015-12-09

**K 7727**

**Immundiagnostik AG**, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: [info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)

## Table of Contents

<b>1. INTENDED USE</b>	<b>14</b>
<b>2. MATERIAL SUPPLIED</b>	<b>14</b>
<b>3. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>14</b>
<b>4. PREPARATION AND STORAGE OF REAGENTS</b>	<b>15</b>
<b>5. PREPARATION AND STORAGE OF SAMPLES</b>	<b>16</b>
<b>6. ASSAY PROCEDURE</b>	<b>16</b>
<i>Principle of the test</i>	16
<i>Test procedure</i>	16
<b>7. RESULTS</b>	<b>18</b>
<b>8. LIMITATIONS</b>	<b>19</b>
<b>9. QUALITY CONTROL</b>	<b>19</b>
<b>10. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
<i>Analytical sensitivity</i>	20
<b>11. PRECAUTIONS</b>	<b>20</b>
<b>12. TECHNICAL HINTS</b>	<b>20</b>
<b>13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>21</b>
<b>14. REFERENCES</b>	<b>21</b>

## 1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of indoleamine 2,3-dioxygenase 1 (IDO1). It is for research use only. Not for use in diagnostic procedures.

## 2. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 7727	PLATE	Holder with precoated strips	12 x 8 Wells
K 7727	STD	Standards lyophilized (0, 0.3, 1.5, 6, 25, 100 ng/ml)	1 x 6 vials
K 7727	CTRL 1 CTRL 2	Controls, lyophilized (see specification for range)	1 x 2 vials
K 7727	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 7727	AB	IDO antibody, lyophilized	2 x 1 vial
K 7727	ASYBUF	Assay buffer, ready to use	1 x 30 ml
K 7727	CONJ	Conjugate (peroxidase-labeled), concentrate	1 x 65 µl
K 7727	CONJBUF	Conjugate stabilizing buffer, ready to use	1 x 13 ml
K 7727	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 7727	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

## 3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Calibrated precision pipets and 10-1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 x g

- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 6)

\* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

#### 4. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used several times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37 °C using a water bath before dilution. The wash buffer concentrate is stable at **2-8 °C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8 °C for one month.**
- The lyophilized **standards (STD)** and **controls (CTRL)** are stable at **2-8 °C** until the expiry date stated on the label. Before use, the standards and controls must be reconstituted with **200 µl of assay buffer (ASYBUF)**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to ensure complete reconstitution. Reconstituted standard and controls can be stored at **-20 °C.**
- Dissolve the the content of one vial of **IDO antibody (AB)** in **5.5 ml of assay buffer (ASYBUF)**. When more than one vial is to be used, combine the contents and mix prior to use. The reconstituted IDO antibody can be stored at **-20°C.**
- Dilute the **peroxidase conjugate (CONJ) 1:201** with conjugate stabilizing buffer (CONJBUF) (**e.g. 60 µl CONJ + 12 ml CONJBUF**, prepare only the required amount). The undiluted POD conjugate is stable at 2-8 °C until the expiry date stated on the label. Diluted POD conjugate can be stored at **2-8 °C for 1 week.**

- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8 °C**.

## 5. PREPARATION AND STORAGE OF SAMPLES

Serum and EDTA plasma are suited for this test system. The samples are analyzed **undiluted**.

Analysis of other sample matrices is possible, e.g. cell culture supernatant or tissue homogenates, but different standard concentrations may be necessary. Please do not hesitate to contact Immundiagnostik AG.

Samples with visible amounts of precipitates should be centrifuged. Store samples frozen at -20 °C.

## 6. ASSAY PROCEDURE

### *Principle of the test*

The assay utilizes the two-site sandwich technique with two selected polyclonal antibodies that bind to human indoleamine 2,3-deoxygenase 1 (IDO1).

Assay standards, controls and samples containing IDO1 are added to the wells of a microplate coated with a high affine polyclonal anti-human IDO1 antibody. During the first incubation step IDO1 in the samples is captured by the immobilized antibody. Then, the biotinylated detection antibody, a polyclonal anti-human IDO1 antibody, is added. In the next step peroxidase labeled streptavidin is added to each well and the following complex is formed: capture antibody – human IDO1 – detection antibody – peroxidase conjugate. Tetramethylbenzidine (TMB) is used as a substrate for the peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the IDO1 concentration of sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from standards. The IDO1 concentration of the samples is determined directly from this curve.

### *Test procedure*

Bring all **reagents and samples to room temperature** (15-30 °C) and mix well.

Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2-8 °C. Strips are stable until the expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

We recommend to carry out the tests in duplicate.

1.	Wash each well <b>5 x</b> with <b>250 µl of diluted wash buffer</b> before use. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
2.	Add 2 x <b>25 µl</b> of <b>standards/ controls/ samples</b> (STD/ CTRL/ SAMPLE) into the respective wells of the microtiter plate.
3.	Add <b>100 µl</b> of <b>assay buffer (ASYBUF)</b> into each well
4.	Cover the plate tightly and incubate for <b>2 hours</b> at <b>room temperature</b> (15-30 °C) on a horizontal <b>shaker</b> (180-240 rpm).
5.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
6.	Add <b>100 µl</b> of <b>IDO antibody (AB)</b> into each well.
7.	Cover the plate tightly and incubate for <b>1 hour</b> at <b>room temperature</b> (15-30 °C) on a horizontal <b>shaker</b> (180-240 rpm).
8.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
9.	Add <b>100 µl</b> of diluted <b>peroxidase conjugate (CONJ)</b> into each well.
10.	Cover the plate tightly and incubate for <b>30 min</b> at <b>room temperature</b> (15-30 °C) on a horizontal <b>shaker</b> (180-240 rpm).

11.	Discard the contents of each well and wash <b>5 times</b> with <b>250 µl of diluted wash buffer</b> . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
12.	Add <b>100 µl</b> of <b>TMB substrate (SUB)</b> into each well.
13.	Incubate for <b>12-18 min</b> at <b>room temperature</b> (15-30 °C) in the dark*.
14.	Add <b>100 µl</b> of <b>stop solution (STOP)</b> into each well, mix thoroughly.
15.	Determine <b>absorption immediately</b> with an ELISA reader at <b>450 nm</b> against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at <b>405 nm</b> against 620 nm (690 nm) as a reference.

\* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

## 7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

### 1. 4-parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

### 2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### 3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.


The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

### Serum and EDTA plasma

The concentrations can be determined directly from the standard curve in ng/ml. **No factor** is needed.

If a sample has been diluted, multiply the obtained result by the dilution factor used.

In the following, an example of a standard curve is given; do not use it for the calculation of your results.

## 8. LIMITATIONS

Samples with an OD higher than the OD of the highest standard should be diluted with assay buffer (ASYBUF) and re-assayed. Please consider this dilution factor when calculating the results.

## 9. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using

appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

## 10. PERFORMANCE CHARACTERISTICS

### *Analytical sensitivity*

The zero-standard was measured 42 times. The detection limit was set as  $B_0 + 2 \text{ SD}$  and estimated to be 0.09 ng/ml.

sample	mean value [OD]	2 x standard deviation (2 x SD)	detection limit [ng/ml]
zero-standard	0.111	0.010	0.09

## 11. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes
- The stop solution consists of sulfuric acid, which is a strong acid. Although diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

## 12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.

- Control Samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

### 14. REFERENCES

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; **8** (3): 289-95.

Cavia-Saiz M, Muñiz Rodríguez P, Llorente Ayala B, García-González M, Coma-Del Corral MJ, García Girón C. The role of plasma IDO activity as a diagnostic marker of patients with colorectal cancer. *Mol Biol Rep.* 2014 Apr; **41**(4):2275-9.

Ciorba MA: Indoleamine 2,3 dioxygenase in intestinal disease. *Curr Opin Gastroenterol.* 2013 Mar; **29**(2):146-52

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology*. 2013 Mar 1; **2**(3):e23428

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol*. 2012 March; **19**(3): 436–442

### Used symbols:

Temperature limitation

Catalogue Number

For research use only

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by