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Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev.* 2014 Mar; **23**(3):461-8

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology.* 2013 Mar 1; **2**(3):e23428

Eussen SJPM, Ueland PM, Vollset SE, Nygård O, Midttun Ø, Sulo G, Tell GS: Kynurenines as predictors of acute coronary events in the Hordaland Health Study. *Int J Cardiol.* 2015 Jun 15; **189**:18-24

Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis.* 2012 Jul; **18**(7):1214-20.

Murray C, Griffin ÉW, O'Loughlin E, Lyons A, Sherwin E, Ahmed S, Stevenson N J, et al. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity.* 2015, **48**: 274-286.

Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Iglund J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of

kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J.* 2013 Sep; **34**(34):2689-96.

Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc.* 2014; **3**:e001094

Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol.* 2013 Sep 30; **168**(2):1435-40

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol.* 2012 March; **19**(3): 436–442

Yan EB, Frugier T, Lim C K, Heng B, Sundaram G, Tan M, Rosenfeld J V, et al. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation.* 2015, **12**(1), 110.

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Manual

IDK® Kynurenine ELISA **high sensitive**



For the in vitro determination of L-kynurenine in rodents (plasma, serum and brain tissue), in cell culture supernatant and CSF

For research use only

Valid from 2016-01-12

REF

K 3728

RUO

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: ++49 6251 70190-0

Fax: ++ 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	18
2. MATERIAL SUPPLIED	18
3. MATERIAL REQUIRED BUT NOT SUPPLIED	19
4. PREPARATION AND STORAGE OF REAGENTS	19
5. PREPARATION AND STORAGE OF SAMPLES	20
6. ASSAY PROCEDURE	21
<i>Principle of the test</i>	21
<i>Sample preparation procedure</i>	22
<i>Test procedure</i>	22
7. RESULTS	24
8. LIMITATIONS	26
9. QUALITY CONTROL	26
<i>Reference Range</i>	26
10. PERFORMANCE CHARACTERISTICS	26
<i>Precision and reproducibility</i>	26
<i>Spiking recovery</i>	27
<i>Dilution recovery</i>	27
<i>Analytical sensitivity</i>	28
<i>Specificity</i>	28
11. PRECAUTIONS	28
12. TECHNICAL HINTS	29
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	29
14. REFERENCES	29

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of L-kynurenine in EDTA plasma, serum and brain tissue of rodents (mouse, rat), and in cell culture supernatant and CSF. For research use only. Not for use in diagnostic procedures.

For human plasma, serum, and urine samples we recommend our *IDK*® Kynurenine ELISA K 7728.

2. MATERIAL SUPPLIED

Catalog No	Label	Kit Components	Quantity
K 3728	PLATE	Holder with precoated strips	12 x 8 wells
K 3728	STD	Standards, ready to use (0, 0.04, 0.16, 0.64, 2.5, 10 µM)	6 x 200 µl
K 3728 K 3728	CTRL 1 CTRL 2	Controls, ready to use (see specification for range)	2 x 200 µl
K 3728	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 3728	AB	L-kynurenine antibody, lyophilized	2 x 1 vial
K 3728	CONJ	Conjugate (peroxidase-labeled), concentrate	1 x 65 µl
K 3728	CONJBUF	Conjugate stabilizing buffer, ready to use	1 x 13 ml
K 3728	REABUF	Reaction buffer, ready to use	1 x 45 ml
K 3728	DER	Derivatization reagent	2 x 12.5 mg
K 3728	DMSO	Dimethylsulfoxide (DMSO)	1 x 2ml
K 3728	ASYBUF	Assay buffer, ready to use	1 x 28 ml
K 3728	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 3728	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipets and 10-1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 14000 x g
- Laboratory balance
- Ultrasonic homogenizer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 6)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

4. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- Dilute the **wash buffer concentrate (WASHBUF)** with ultra pure water **1:10** before use (**100 ml WASHBUF + 900 ml ultra pure water**), mix well. Crystals may occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or at 37°C using a water bath before dilution. The wash buffer concentrate is stable at **2-8 °C** until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at **2-8 °C for one month.**
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20 °C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use. They can be re-frozen up to 3 times.

- **DMSO** could crystallize at 4 °C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER)** (12.5 mg) in **750 µl DMSO**. Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use**. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. The ELISA kit can be separated into two performances by providing two DER vials. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.
- Dissolve the content of one vial of **L-kynurenine antibody (AB)** in **3 ml of diluted wash buffer**. When more than one vial is to be used, combine the contents and mix prior to use. The ELISA kit can be separated into two performances by providing two AB vials. Diluted L-kynurenine antibody can be stored at **-20 °C for one month**.
- Dilute the **peroxidase conjugate (CONJ) 1:201** with conjugate stabilizing buffer (CONJBUF) (**e.g. 60 µl CONJ + 12 ml CONJBUF**, prepare only the required amount). The undiluted POD conjugate is stable at 2-8 °C until the expiry date stated on the label. Diluted POD conjugate can be stored at **2-8 °C for 1 week**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8 °C**.

5. PREPARATION AND STORAGE OF SAMPLES

EDTA plasma, serum, cell culture supernatant, CSF

- EDTA-Plasma and serum of rodents (mouse, rat), cell culture supernatant and CSF (cerebrospinal fluid) is suited for this test system. The samples are stable for 3 days at 2-8 °C. For longer storage samples must be frozen at -20 °C or -80 °C.
- Samples with visible amounts of precipitates should be centrifuged.
- Samples are used **undiluted**.

Brain tissue

- Brain tissue of rodents (mouse, rat) can be tested. Place the tissue sample in a microcentrifuge tube and determine the exact wet mass. Add 40 µl of reaction buffer (REABUF) per milligram of tissue to the tube, e.g.

25 mg tissue + 1000 µl REABUF. Homogenize the sample by sonication and centrifuge at 14,000 x g for 10 min at 4 °C. Transfer the supernatant to a clean tube. 100 µl of the supernatant are used in the test (see sample preparation procedure). Store samples and supernatants at -20 °C or -80 °C.

For sample preparation, a derivatization reagent (DER) for derivatization of L-kynurenine is added to all samples (details are given in the sample preparation procedure).

6. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for L-kynurenine derivatization. Afterwards, the treated samples and a polyclonal L-kynurenine-antiserum are incubated in wells of a microtiter plate coated with L-kynurenine-derivative (tracer).

During the incubation period the target L-kynurenine in the sample competes with the tracer, immobilized on the wall of the microtiter wells, for the binding of the polyclonal antibodies. The L-kynurenine in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-bound antibody is inverse proportional to the kynurenine concentration in the sample.

During the second incubation step, a peroxidase-conjugated antibody is added to each microtiter well to detect the anti-kynurenine antibodies. After washing away the unbound components, tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the L-kynurenine concentration in the sample; this means, high L-kynurenine concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. L-kynurenine present in the patient samples is determined directly from this curve.

Sample preparation procedure

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials).

The reagents provided with this kit are sufficient for up to 48 derivatizations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1. Bring all reagents and samples to **room temperature** (15-30 °C) and mix well.

2. *Standards/controls, plasma, serum, cell culture supernatant, CSF:*
Add **10 µl** of standards (STD), controls (CTRL), samples (SAMPLE) in the corresponding vials.
Add **100 µl** of **reaction buffer (REABUF)** into each vial (STD, CTRL, SAMPLE).

Brain tissue:

Add **100 µl** of supernatant after centrifugation into the corresponding vial (SAMPLE).

3. Add **25 µl** of freshly prepared **derivatization reagent (DER)** into each vial (STD, CTRL, SAMPLE) and **mix thoroughly** by repeated inversion or several seconds on a vortex mixer. Incubate for **45 min at room temperature** (15-30 °C) on a horizontal **shaker** (180-240 rpm).

4. Afterwards add **250 µl of assay buffer (ASYBUF)** into each vial, **mix thoroughly** by repeated inversion or on a vortex mixer and incubate for **5 min at room temperature** (15-30 °C) on a horizontal **shaker** (180-240 rpm).

2 x 50 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5. Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) in duplicate on a **protocol sheet**.

6. Take as many microtiter strips (PLATE) as needed from kit. Store unused strips covered at 2-8 °C. Strips are stable until the expiry date stated on the label.
7. Wash each well 5 x with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
8. For the analysis in duplicate, take 2 x 50 µl of the derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE) out of the vials and add into the respective wells of the microtiter plate.
9. Add 50 µl of L-kynurenine antibody (AB) into each well. Cover the plate tightly.
10. Incubate overnight (15-20 hours) at 2-8 °C .
11. Discard the contents of each well. Wash each well 5 x with 250 µl of diluted wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
12. Add 100 µl of diluted peroxidase conjugate (CONJ) into each well.
13. Cover the plate tightly and incubate for 1 hour at room temperature (15-30°C) on a horizontal shaker (180-240 rpm).
14. Discard the contents of each well. Wash each well 5 x with 250 µl of diluted wash buffer . After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
15. Add 100 µl of TMB substrate (SUB) into each well.
16. Incubate for 8-12 min at room temperature (15-30 °C) in the dark*.
17. Add 100 µl of stop solution (STOP) into each well, mix thoroughly.

18. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm (690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4-parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

EDTA plasma, serum, CSF

The concentrations can be determined directly from the standard curve in $\mu\text{mol/l}$. No factor is needed.

Cell culture supernatant

Cell culture media show matrix effects when very low levels have to be determined.


We suggest you proceed as follows: High L-kynurenine levels from 1 to 10 $\mu\text{mol/l}$ in the sample can be determined with the standard curve provided with the kit. For samples with lower L-kynurenine levels we recommend preparing an own standard curve by diluting the highest standard (STD6) with the cell culture media used in your laboratory.

The concentrations can be determined directly from the standard curve. No factor is needed.

Brain tissue

The obtained L-kynurenin levels of brain tissue samples have to be multiplied by the **factor of 4**.

In the following, an example of a standard curve is given; do not use it for the calculation of your results.

8. LIMITATIONS

Samples with an OD lower than the OD of the highest standard should be diluted with reaction buffer (REABUF) and re-assayed. Please consider this dilution factor when calculating the results.

9. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference Range

CSF: Yan, E.B. et al., 2015. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*, 12(1):110.

Plasma (mice): Murray, C. et al., 2015. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity*, 48:274-286.

We recommend each laboratory to establish its own reference range.

10. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 10)

Sample	L-kynurenine [$\mu\text{mol/l}$]	CV [%]
1	1.04	8.1
2	2.94	4.9

Inter-Assay (n = 10)

Sample	L-kynurenine [$\mu\text{mol/l}$]	CV [%]
1	1.07	9.3
2	2.83	11.6

Spiking recovery

Two samples were spiked with different L-kynurenine concentrations and measured in this assay. The mean recovery rate was 98.7 % (n = 2).

Spike [$\mu\text{mol/l}$]	L-kynurenine expected [$\mu\text{mol/l}$]	L-kynurenine measured [$\mu\text{mol/l}$]	Recovery [%]
sample 1		0.98	
1.5	2.48	2.25	90.1
3.0	3.98	3.95	99.2
sample 2		1.13	
1.5	2.63	2.51	95.4
3.0	4.13	4.21	102.0

Dilution recovery

Two spiked samples were diluted with reaction buffer and measured in this assay. The mean recovery rate was 97.4 % (n = 2).

Dilution	L-kynurenine expected [$\mu\text{mol/l}$]	L-kynurenine measured [$\mu\text{mol/l}$]	Recovery [%]
sample 1		2.25	
1:2	1.13	1.19	105.7
1:4	0.56	0.55	97.2
sample 2		3.95	
1:2	1.97	1.85	94.0
1:4	0.99	0.92	92.8

Analytical sensitivity

The zero-standard was measured 48 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be 0.025 [$\mu\text{mol/l}$]. Considering the dilution factor, the detection limit for brain tissue samples is calculated to be 0.1 [$\mu\text{mol/l}$].

Zero standard mean value [OD]	2 x standard deviation (SD)	detection limit [$\mu\text{mol/l}$]	detection limit brain tissue [$\mu\text{mol/l}$]
1.941	0.184	0.025	0.1

Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to L-kynurenine. The specificity is calculated in percent in relation to the L-kynurenine binding activity.

3-HK (3-hydroxy-DL-kynurenine)	< 0.5 %
L-tryptophan	< 0.08 %

11. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes
- The stop solution consists of sulfuric acid, which is a strong acid. Although diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

12. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control Samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

14. REFERENCES

Cavia-Saiz M, Muñoz Rodríguez P, Llorente Ayala B, García-González M, Coma-Del Corral MJ, García Girón C. The role of plasma IDO activity as a diagnostic marker of patients with colorectal cancer. *Mol Biol Rep*. 2014 Apr; **41**(4):2275-9.

Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev*. 2014 Mar; **23**(3):461-8

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology*. 2013 Mar 1; **2**(3):e23428

Eussen S J PM, Ueland P M, Vollset S E, Nygård O, Midttun Ø, Sulo G, Tell GS: Kynurenines as predictors of acute coronary events in the Hordaland Health Study. *Int J Cardiol*. 2015 Jun 15; **189**:18-24

Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis*. 2012 Jul; **18**(7):1214-20.

Murray C, Griffin ÉW, O'Loughlin E, Lyons A, Sherwin E, Ahmed S, Stevenson N J, et al. Interdependent and independent roles of type I interferons and IL-6 in innate immune, neuroinflammatory and sickness behaviour responses to systemic poly I:C. *Brain, Behavior, and Immunity*. 2015, **48**: 274-286.

Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Igland J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J*. 2013 Sep; **34**(34):2689-96.

Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc*. 2014; **3**:e001094

Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol*. 2013 Sep 30; **168**(2):1435-40

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol*. 2012 March; **19**(3): 436–442

Yan EB, Frugier T, Lim CK, Heng B, Sundaram G, Tan M, Rosenfeld JV, et al. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*. 2015, **12**(1), 110.

Used symbols:

Temperature limitation

Catalogue Number

For research use only

To be used with

Manufacturer

Contains sufficient for <n> tests

Lot number

Use by

Attention