ClinCal® - Calibrator

Dried Blood Spot Calibrator Set / Trockenblut-Kalibrator-Set (Filterspots)

FOR HOMOCYSTEINE / FÜR HOMOCYSTEIN

Intended use:

ClinCal® dried blood spot calibrators are used for calibration of the analysis-system. These calibrators are based on human whole blood and are available as a calibrator set with 4 levels.

Pretreatment:

After having stamped a defined punching out of the blood spot, it has to be prepared like patient samples in one series of analyses. Please note that the blood volume that was used for creating one spot on the filter paper is exactly **100 µl**. The diameter of one filter spot is **18 mm**. If you use punchings of smaller diameters, you have to calculate the corresponding blood volume of it (**area ratio!**).

Storage and stability:

Originally closed and stored at 2 - 8 °C, the dried blood spot calibrator set is stable for 36 months, but not longer than the expiration date printed on the label.

Before opening the cooled calibrator pack, you have to temper it to room temperature, to protect against condensed water. The pack includes a pouch with desiccant (orange gel), which serves as indicator for moisture simultaneously. In the dry status the gel is coloured orange. If the capacity for water absorption is expended, the desiccant will change to blue. Then the material must be regenerated by heating for at least 5 hours at 70 °C in a cabinet dryer. After cooling down, it should be enclosed to the filter cards in the alumina bag.

Opened calibrator packs can be stored until the expiration date, printed on the label, if they are stored dry and at 2 - 8 °C.

Notes:

The mean values of the analyte are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each filter test card of one lot. The variation (CV) of the volume per blood spot is < 1 %.

Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. The determination of the analyte was performed by LC-MS/MS.

Pack size:

ClinCal® Dried Blood Spot Calibrator Set Level 0 - 3
4 x 1 x 3 spots, **order no.: MS2113**

Precautions:

The human whole blood which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1-, HCV-RNA and HBV-DNA (NAT). Nevertheless, the calibrators should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH