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IDK[®] Tryptophan ELISA

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15. LITERATUR

Allgemeine Literatur

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; **8** (3): 289-95.

Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev.* 2014 Mar;**23**(3):461-8

Ciorba MA: Indoleamine 2,3 dioxygenase in intestinal disease. *Curr Opin Gastroenterol.* 2013 Mar;**29**(2):146-52

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology.* 2013 Mar 1;**2**(3):e23428

Dolina S, Margalit D, Malitsky S, Rabinkov A: Attention-deficit hyperactivity disorder (ADHD) as a pyridoxine-dependent condition: urinary diagnostic biomarkers. *Med Hypotheses.* 2014 Jan;**82**(1):111-6.

Grozdics E, Berta L, Bajnok A, Veres G, Ilisz I, Klivényi P, Rigó J Jr, Vécsei L, Tulassay T, Toldi G: B7 costimulation and intracellular indoleamine-2,3-dioxygenase (IDO) expression in peripheral blood of healthy pregnant and non-pregnant women. *BMC Pregnancy Childbirth.* 2014 Sep 4;**14**:306

Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis.* 2012 Jul;**18**(7):1214-20.

Hyangin Kim, Lucy Chen, Grewo Lim, Backil Sung, Shuxing Wang, Michael F. McCabe, Gabriel Rusanescu, Liling Yang, Yinghong Tian, Jianren Mao: Brain indoleamine 2,3-dioxygenase contributes to the comorbidity of pain and depression. *J Clin Invest.* 2012 August 1; **122**(8): 2940–2954.

Miura H, Ozaki N, Sawada M, Isobe K, Ohta T, Nagatsu T: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress.* 2008 May;**11**(3):198-209.

Myint AM, Bondy B, Baghai TC, Eser D, Nothdurfter C, Schüle C, Zill P, Müller N, Rupprecht R, Schwarz MJ: Tryptophan metabolism and immunogenetics in

major depression: a role for interferon- γ gene. *Brain Behav Immun*. 2013 Jul;**31**:128-33

Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Igland J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J*. 2013 Sep;**34**(34):2689-96.

Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc*. 2014;**3**:e001094

Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol*. 2013 Sep 30;**168**(2):1435-40

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol*. 2012 March; **19**(3): 436–442

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *AM J Med Genet B Neuropsychiatr Genet*. 2009 Mar 5; **150B** (2): 202-8.

Yan EB, Frugier T, Lim C K, Heng B, Sundaram G, Tan M, Rosenfeld J V, et al.: Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*. 2015,**12**(1), 110.

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
Manual

IDK® Tryptophan ELISA

For the in vitro determination of L-tryptophan in human EDTA plasma, serum and urine

Valid from 2016-08-22

REF K 7730

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	21
2. INTRODUCTION	21
3. MATERIAL SUPPLIED	22
4. MATERIAL REQUIRED BUT NOT SUPPLIED	23
5. PREPARATION AND STORAGE OF REAGENTS	23
6. PREPARATION AND STORAGE OF SAMPLES	24
7. ASSAY PROCEDURE	25
<i>Principle of the test</i>	25
<i>Sample preparation procedure</i>	25
<i>Test procedure</i>	26
8. RESULTS	27
9. LIMITATIONS	28
10. QUALITY CONTROL	29
<i>Reference Range</i>	29
11. PERFORMANCE CHARACTERISTICS	30
<i>Precision and reproducibility</i>	30
<i>Spiking recovery</i>	30
<i>Dilution recovery</i>	31
<i>Analytical sensitivity</i>	31
<i>Specificity</i>	31
<i>Correlation with HPLC-MS</i>	32
12. PRECAUTIONS	32
13. TECHNICAL HINTS	33
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	33
15. REFERENCES	34
<i>General Literature</i>	34
<i>Publications using Immundiagnostik IDK® Tryptophan ELISA</i>	35

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of L-tryptophan in human EDTA plasma, serum and urine. For *in vitro* diagnostic use only.

For rodent specimens (mouse, rat) and for cell culture supernatant and CSF we recommend our IDK® Tryptophan high sensitive ELISA K 3730.

2. INTRODUCTION

Tryptophan is an essential amino acid which cannot be synthesized by humans, and therefore must be supplied in diet. Tryptophan plays an important role in the synthesis of serotonin, as it is the only source for synthesizing serotonin (see figure 1).

Serotonin is released from the appendices of a widespread transmitter system with global modulating effects on different areas of life. The activity of this system depends not least on the actual tryptophan level.


Figure 1: Synthesis of serotonin and melatonin *via* tryptophan (left), and shifting of the metabolism of tryptophan to the kynurenine-pathway *via* the enzyme indoleamine 2,3-dioxygenase (right).

When tryptophan in plasma decreases, it is not available for serotonin synthesis and the following melatonin production. As a result, fluctuations of mood, aggression, sleep and eating disorder, or depression can occur. The symptoms of these diseases can be eased by administration of tryptophan. However, such a treatment requires the knowledge of the tryptophan level, which can be determined with this test.

3. MATERIAL SUPPLIED

Catalog No	Label	Kit Components	Quantity
K 7730	PLATE	Holder with precoated strips	12 x 8 wells
K 7730	STD	Standards, lyophilized (see specification for concentrations)	1 x 6 vials
K 7730	CTRL 1	Control, lyophilized (see specification for range)	1 vial
K 7730	CTRL 2	Control, lyophilized (see specification for range)	1 vial
K 7730	WASHBUF	ELISA wash buffer concentrate, 10x	2 x 100 ml
K 7730	AB	L-tryptophan antibody, lyophilized	4 x 1 vial
K 7730	CONJ	Conjugate (peroxidase-labeled), concentrate	1 x 65 µl
K 7730	CONJBUF	Conjugate stabilizing buffer, ready to use	1 x 13 ml
K 7730	REABUF	Reaction buffer, ready to use	2 x 28 ml
K 7730	DER	Derivatization reagent	4 x 25 mg
K 7730	DMSO	Dimethylsulfoxide (DMSO)	1 x 7 ml
K 7730	SUB	TMB substrate (tetramethylbenzidine), ready to use	1 x 15 ml
K 7730	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipets and 10-1000 µl tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 x g
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month**.
- The **lyophilized standards (STD) and controls (CTRL)** are stable at **2-8 °C** until the expiry date stated on the label. Before use, reconstitute the standards and controls with **200 µl of reaction buffer (REABUF)**. Mix thoroughly and allow to dissolve for **15 min** on a horizontal shaker. **Reconstituted standards and controls can be stored at -20 °C.**

- **DMSO** crystallizes at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (25 mg) in 1.5 ml DMSO**. Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use**. When more than one vial is to be used, combine the contents and mix prior to use. Discard any rest of the reagent after use. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.
- Dissolve the content of one vial of **L-tryptophan antibody (AB) in 3 ml of wash buffer**. When more than one vial is to be used, combine the contents and mix prior to use. Diluted L-tryptophan antibody can be stored at **-20 °C for 1 month**.
- Dilute the **peroxidase conjugate (CONJ) 1:201** with **conjugate stabilizing buffer (CONJBUF)** (e.g. 60 µl CONJ + 12 ml CONJBUF, prepare only the required amount). The undiluted peroxidase conjugate is stable at 2-8 °C until the expiry date stated on the label. Diluted peroxidase conjugate can be stored at **2-8 °C for 1 week**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8 °C**.

6. PREPARATION AND STORAGE OF SAMPLES

EDTA plasma, serum and urine

- Serum, EDTA-Plasma and urine are suited for this test system. The samples are stable for 3 days at 2-8°C and at room temperature. For longer storage keep samples frozen at -20 °C. We recommend acidifying the urine samples.
- Samples are used **undiluted**.
Samples from patients with tryptophan supplementation (e.g. in depletion studies) probably require dilution, we recommend diluting these samples 1:2 with REABUF (e.g. 25 µl sample + 25 µl REABUF).
- For sample preparation, a derivatization reagent (DER) for derivatization of tryptophan is added (details are given in the sample preparation procedure).

7. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for L-tryptophan derivatization. Afterwards, the treated samples and a polyclonal L-tryptophan antiserum are incubated in the wells of a microtiter plate coated with L-tryptophan derivative (tracer).

During the incubation period, the target L-tryptophan in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the polyclonal antibodies. The L-tryptophan in the sample displaces the antibodies out of the binding to the tracer. Therefore the concentration of the tracer-bound antibody is inverse proportional to the L-tryptophan concentration in the sample.

During the second incubation step a peroxidase-conjugated antibody is added to each microtiter well to detect the L-tryptophan antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as a peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the L-tryptophan concentration in the sample; this means, high L-tryptophan concentration in the sample reduces the concentration of tracer-bound antibodies and lowers the photometric signal. A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. L-tryptophan present in the patient samples is determined directly from this curve.

Sample preparation procedure

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials).

We recommend carrying out 48 derivatizations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1.	Bring all reagents and samples to room temperature (15-30 °C) and mix well.
----	--

2.	Add 25 µl of standards (STD), controls (CTRL) or samples (SAMPLE) in the corresponding vials.
3.	Add 500 µl of reaction buffer (REABUF) into each vial (STD, CTRL, SAMPLE).
4.	Add 50 µl of freshly prepared derivatization reagent (DER) into each vial (STD, CTRL, SAMPLE) and mix thoroughly by repeated inversion or several seconds on a vortex mixer. Incubate for 45 min at room temperature (15-30°C) on a horizontal shaker .

2 x 25 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5.	Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) in duplicate on a protocol sheet .
6.	Take as many microtiter plate strips (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label. Please note: Do not wash the plate .
7.	For the analysis in duplicate, take 2 x 25 µl of the derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE) out of the vials and add into the respective wells of the microtiter plate.
8.	Add 100 µl of L-tryptophan antibody (AB) into each well. Cover the plate tightly.
9.	Incubate for 2 hours at room temperature (15-30°C) on a horizontal shaker .
10.	Discard the content of each well and wash 5 times with 250 µl of wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
11.	Add 100 µl of diluted peroxidase conjugate into each well.

12.	Cover the plate and incubate for 1 hour at room temperature (15-30 °C) on a horizontal shaker .
13.	Discard the content of each well and wash 5 times with 250 µl of wash buffer . After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
14.	Add 100 µl of TMB substrate (SUB) into each well.
15.	Incubate for 10-15 min at room temperature (15-30 °C) in the dark*.
16.	Add 100 µl of ELISA stop solution (STOP) into each well, mix thoroughly.
17.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4-parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.


The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.

EDTA plasma, serum, urine

The concentrations can be determined directly from the standard curve in $\mu\text{mol/l}$. **No factor** is required.

The results from diluted samples must be multiplied by this dilution factor.

In the following, an example of a standard curve is given; do not use it for the calculation of your results.

9. LIMITATIONS

Samples with concentrations above the measurement range (see definition below) must be diluted with REABUF and re-assayed. Please consider this dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

highest concentration of the standard curve × sample dilution factor to be used

The lower limit of the measurement range can be calculated as:

LoB × sample dilution factor to be used

10. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference Range

Serum

Based on internal studies with serum samples of apparently healthy persons (n = 70) a mean value of 63.5 µmol/l was estimated. The standard deviation was 10.1 µmol/l.

Mean value ± 2x standard deviation: **63.5 ± 20.2 µmol/l**

Normal range: **43.3 – 83.7 µmol/l**

We recommend each laboratory to establish its own reference range.

Urine

Based on internal studies with urine samples of apparently healthy persons (n = 42) the following reference range was determined:

Median: **54.6 µmol/g creatinine**

For 95 % of this collective an L-tryptophan concentration in the range of 29 - 165 µmol/g creatinine was obtained.

We recommend each laboratory to establish its own reference range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-assay (n = 14)

sample	L-tryptophan [$\mu\text{mol/l}$]	CV [%]
1	51.4	4.3
2	105.7	6.7

Inter-assay (n = 7)

sample	L-tryptophan [$\mu\text{mol/l}$]	CV [%]
1	63.7	8.4
2	60.6	9.1

Spiking recovery

3 serum sample were spiked with different tryptophan concentrations and measured in this assay. The mean recovery rate for all concentrations was 97.2 % (n = 2).

spike [$\mu\text{mol/l}$]	L-tryptophan expected [$\mu\text{mol/l}$]	L-tryptophan measured [$\mu\text{mol/l}$]	recovery [%]
		66.8	
50	116.8	108.0	92.5
100	166.8	151.1	90.6
		64.1	
50	114.1	122.1	107.0
100	164.1	184.9	112.7
		69.9	
50	119.9	112.8	94.2
100	169.9	146.9	86.5

Dilution recovery

2 serum samples were diluted and measured in this assay. The mean recovery was 96.8 % (n = 2).

dilution	L-tryptophan expected [μmol/l]	L-tryptophan measured [μmol/l]	recovery [%]
		91.5	
1:2	45.8	46.7	102.0
1:3	30.5	27.5	90.1
1:4	22.9	22.7	99.4
1:5	18.3	17.7	96.6
		89.6	
1:2	44.8	41.2	91.9
1:3	29.9	25.2	84.5
1:4	22.4	21.9	97.6
1:5	17.9	21.1	117.5

Analytical sensitivity

Limit of blank, LoB 4.9 μmol/l

Limit of detection, LoD 8.0 μmol/l

Limit of quantitation, LoQ 10.0 μmol/l

The evaluation was performed according to the CLSI guideline EP-17-A2. The specified accuracy goal for the LoQ was 15 % CV.

Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to tryptophan. The specificity is calculated in percent in relation to the tryptophan binding activity.


5-HTP (5-Hydroxy-Tryptophan) < 0.5 %

L-phenylalanine < 0.1 %

L-tyrosine < 0.1 %

Correlation with HPLC-MS

10 samples were measured with this ELISA and HPLC-MS. The correlation was $r = 0.98$.

12. PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes
- The stop solution consists of sulfuric acid, which is a strong acid. Although diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control Samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- *IDK®* is a trademark of Immundiagnostik AG.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

15. REFERENCES

General Literature

Brandacher G, Hoeller E, Fuchs D, Weiss HG. Chronic immune activation underlies morbid obesity: is IDO a key player. *Curr Drug Metab.* 2007; **8** (3): 289-95.

Chuang SC, Fanidi A, Ueland PM et al: Circulating biomarkers of tryptophan and the kynurenine pathway and lung cancer risk. *Cancer Epidemiol Biomarkers Prev.* 2014 Mar;**23**(3):461-8

Ciorba MA: Indoleamine 2,3 dioxygenase in intestinal disease. *Curr Opin Gastroenterol.* 2013 Mar;**29**(2):146-52

Creelan BC, Antonia S, Bepler G, Garrett TJ, Simon GR, Soliman HH: Indoleamine 2,3-dioxygenase activity and clinical outcome following induction chemotherapy and concurrent chemoradiation in Stage III non-small cell lung cancer. *Oncoimmunology.* 2013 Mar 1;**2**(3):e23428

Dolina S, Margalit D, Malitsky S, Rabinkov A: Attention-deficit hyperactivity disorder (ADHD) as a pyridoxine-dependent condition: urinary diagnostic biomarkers. *Med Hypotheses.* 2014 Jan;**82**(1):111-6.

Grozdics E, Berta L, Bajnok A, Veres G, Ilisz I, Klivényi P, Rigó J Jr, Vécsei L, Tulassay T, Toldi G: B7 costimulation and intracellular indoleamine-2,3-dioxygenase (IDO) expression in peripheral blood of healthy pregnant and non-pregnant women. *BMC Pregnancy Childbirth.* 2014 Sep 4;**14**:306

Gupta NK, Thaker AI, Kanuri N, Riehl TE, Rowley CW, Stenson WF, Ciorba MA: Serum analysis of tryptophan catabolism pathway: correlation with Crohn's disease activity. *Inflamm Bowel Dis.* 2012 Jul;**18**(7):1214-20.

Hyangin Kim, Lucy Chen, Grewo Lim, Backil Sung, Shuxing Wang, Michael F. McCabe, Gabriel Rusanescu, Liling Yang, Yinghong Tian, Jianren Mao: Brain indoleamine 2,3-dioxygenase contributes to the comorbidity of pain and depression. *J Clin Invest.* 2012 August 1; **122**(8): 2940–2954.

Miura H, Ozaki N, Sawada M, Isobe K, Ohta T, Nagatsu T: A link between stress and depression: Shifts in the balance between the kynurenine and serotonin pathways of tryptophan metabolism and the etiology and pathophysiology of depression; *Stress.* 2008 May;**11**(3):198-209.

Myint AM, Bondy B, Baghai TC, Eser D, Nothdurfter C, Schüle C, Zill P, Müller N, Rupprecht R, Schwarz MJ: Tryptophan metabolism and immunogenetics in

major depression: a role for interferon- γ gene. *Brain Behav Immun*. 2013 Jul;**31**:128-33

Pedersen ER, Svingen GF, Schartum-Hansen H, Ueland PM, Ebbing M, Nordrehaug JE, Igland J, Seifert R, Nilsen RM, Nygård O: Urinary excretion of kynurenine and tryptophan, cardiovascular events, and mortality after elective coronary angiography. *Eur Heart J*. 2013 Sep;**34**(34):2689-96.

Ristagno G, Latini R, Vaahersalo J, Masson S, Kurola J, Varpula T, Lucchetti J, Fracasso C, Guiso G, Montanelli A, Barlera S, Gobbi M, Tiainen M, Pettilä V, Skrifvars MB; FINNRESUSCI Investigators: Early activation of the kynurenine pathway predicts early death and long-term outcome in patients resuscitated from out-of-hospital cardiac arrest. *J Am Heart Assoc*. 2014;**3**:e001094

Sulo G, Vollset SE, Nygård O, Midttun Ø, Ueland PM, Eussen SJ, Pedersen ER, Tell GS: Neopterin and kynurenine-tryptophan ratio as predictors of coronary events in older adults, the Hordaland Health Study. *Int J Cardiol*. 2013 Sep 30;**168**(2):1435-40

Suzuki Y, Suda T, Asada K, Miwa S, Suzuki M, Fujie M, Furuhashi K, Nakamura Y, Inui N, Shirai T, Hayakawa H, Nakamura H, Chida K: Serum Indoleamine 2,3-Dioxygenase Activity Predicts Prognosis of Pulmonary Tuberculosis. *Clin Vaccine Immunol*. 2012 March; **19**(3): 436–442

Wilson ST, Stanley B, Brent DA, Oquendo MA, Huang YY, Mann JJ: The tryptophan hydrolase-1 A218C polymorphism is associated with diagnosis, but not suicidal behavior, in borderline personality disorder. *AM J Med Genet B Neuropsychiatr Genet*. 2009 Mar 5; **150B** (2): 202-8.

Yan EB, Frugier T, Lim C K, Heng B, Sundaram G, Tan M, Rosenfeld J V, et al. Activation of the kynurenine pathway and increased production of the excitotoxin quinolinic acid following traumatic brain injury in humans. *Journal of Neuroinflammation*. 2015,**12**(1), 110.

Publications using Immundiagnostik IDK® Tryptophan ELISA

Hildebrand P, Königschulte W, Gaber TJ, Bubenzer-Busch S, Helmbold K, Biskup CS, Langen KJ, Fink GR, Zepf FD: Effects of dietary tryptophan and phenylalanine-tyrosine depletion on phasic alertness in healthy adults - A pilot study. *Food & nutrition research*. 2015;**59**:26407. doi:10.3402/fnr.v59.26407.

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		