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Tyrosin ELISA

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Manual

Tyrosine ELISA



*For the in vitro determination of tyrosine in human EDTA plasma
and serum*

For research use only

Valid from 2016-06-16

K 7715

Immundiagnostik AG, Stubenwald-Allee 8a, 64625 Bensheim, Germany

Tel.: +49 6251 70190-0

Fax: + 49 6251 849430

e.mail: info@immundiagnostik.com

www.immundiagnostik.com

Table of Contents

1. INTENDED USE	17
2. MATERIAL SUPPLIED	17
3. MATERIAL REQUIRED BUT NOT SUPPLIED	17
4. PREPARATION AND STORAGE OF REAGENTS	18
5. PREPARATION AND STORAGE OF SAMPLES	19
6. ASSAY PROCEDURE	19
<i>Principle of the test</i>	19
<i>Sample preparation procedure</i>	20
<i>Test procedure</i>	21
7. RESULTS	22
8. LIMITATIONS	23
9. QUALITY CONTROL	24
<i>Reference Range</i>	24
10. PERFORMANCE CHARACTERISTICS	24
<i>Precision and reproducibility</i>	24
<i>Spiking recovery</i>	25
<i>Dilution recovery</i>	25
<i>Analytical sensitivity</i>	25
<i>Specificity</i>	26
11. PRECAUTIONS	26
12. TECHNICAL HINTS	26
13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	27

1. INTENDED USE

This Immundiagnostik assay is intended for the quantitative determination of tyrosine in human EDTA plasma and serum. For research use only. Not for use in diagnostic procedures.

2. MATERIAL SUPPLIED

Catalog No	Label	Kit Components	Quantity
K 7015	PLATE	Holder with precoated strips	12 x 8 wells
K 7015	STD	Standards, ready to use (0, 6, 20, 60, 200, 600 μ M)	6 x 200 μ l
K 7015 K 7015	CTRL 1 CTRL 2	Controls, ready to use (see specification for range)	2 x 200 μ l
K 7015	WASHBUF	Wash buffer concentrate, 10x	2 x 100 ml
K 7015	AB	Tyrosine antibody, lyophilized	2 x 1 vial
K 7015	2.AB	Peroxidase labeled antibody (concentrate, 200x)	1 x 90 μ l
K 7015	2.ABDIL	Conjugate stabilizing buffer, ready to use	1 x 16 ml
K 701	REABUF	Reaction buffer, ready to use	1 x 100 ml
K 7015	DER	Derivatization reagent	2 x 13.3 mg
K 7015	DMSO	Dimethylsulfoxide (DMSO)	1 x 2 ml
K 7015	ASYBUF	Assay buffer, ready to use	2 x 28 ml
K 7015	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 7015	STOP	ELISA stop solution, ready to use	1 x 15 ml

For reorders of single components, use the catalogue number followed by the label as product number.

3. MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water*
- Calibrated precision pipets and 10-1000 μ l tips

- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Vortex
- Centrifuge, 3000 g
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter6)

* Immundiagnostik AG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥18.2 MΩ cm).

4. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 2 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:** The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 ml WASHBUF + 900 ml ultra pure water), mix well. Crystals could occur due to high salt concentration in the stock solution. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8 °C** until the expiry date stated on the label. **Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2-8 °C for 1 month.**
- Store **Standards (STD) and controls (CTRL1, CTRL2)** frozen at **-20°C**, thaw before use in the test and mix well. Re-freeze standards and controls immediately after use.
- **DMSO** could crystallize at 4°C. Dissolve the crystals at room temperature or in a water bath.
- Dissolve the content of one vial of **derivatization reagent (DER) (13.3 mg) in 0.8 ml DMSO.** Put the vial on a horizontal shaker for 5 min. DER must be **prepared immediately before use.** When more than one vial is to be used,

combine the contents and mix prior to use. Discard any rest of the reagent after use. **Please note:** DMSO attacks all plastics but not polypropylene products and laboratory glass.

- Dissolve the content of one vial of **tyrosine antibody (AB)** in **5.5 ml of diluted wash buffer**. When more than one vial is to be used, combine the contents and mix prior to use. Diluted anti-tyrosine antibody can be stored at **2-8°C for 1 month**.
- Dilute the **peroxidase labeled antibody (2.AB) 1:201** with **conjugate stabilizing buffer (2.ABDIL)** (e.g. 60 µl 2.AB + 12 ml 2.ABDIL, prepare only the required amount). The undiluted peroxidase labeled antibody is stable at 2-8°C until the expiry date stated on the label. Diluted peroxidase labeled antibody can be stored at **2-8°C for 1 week**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

5. PREPARATION AND STORAGE OF SAMPLES

Serum and EDTA plasma

- Serum and EDTA-Plasma samples are stable for up to 6 hours at 2-8°C. For longer storage samples must be frozen at -20°C.
- The EDTA plasma and serum samples are diluted for derivatization (see sample preparation procedure).
- For sample preparation, a derivatization reagent (DER) for derivatization of tyrosine is added (details are given in the sample preparation procedure).

6. ASSAY PROCEDURE

Principle of the test

This assay is based on the method of competitive enzyme linked immunoassays. The sample preparation includes the addition of a derivatization reagent for tyrosine derivatization.

Afterwards, the treated samples and a polyclonal tyrosine antiserum are incubated in wells of a microtiter plate coated with tyrosine-derivative (tracer). During the incubation period, the target tyrosine in the sample competes with the tracer immobilized on the wall of the microtiter wells for the binding of the

polyclonal antibodies. The tyrosine in the sample displaces the antibodies out of the binding to the tracer. Therefore, the concentration of the tracer-bound antibody is inverse proportional to the tyrosine concentration in the sample.

During the second incubation step, a peroxidase conjugate is added to each microtiter well to detect the anti-tyrosine antibodies. After washing away the unbound components tetramethylbenzidine (TMB) is added as peroxidase substrate. Finally, the enzymatic reaction is terminated by an acidic stop solution. The color changes from blue to yellow, and the absorbance is measured in a photometer at 450 nm. The intensity of the yellow color is inverse proportional to the tyrosine concentration in the sample; this means, high tyrosine concentration in the sample reduces the concentration of tracer-bound antibody and lowers the photometric signal.

A dose response curve of absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Tyrosine present in the patient samples is determined directly from this curve.

Sample preparation procedure

Dilute standards (STD), controls (CTRL) and samples (SAMPLE) by factor **1:41** as follows:

25 µl STD/CTRL/SAMPLE + 1000 µl REABUF (reaction buffer)

Please note: Samples from patients with tyrosine supplementation (e.g. in depletion studies) probably require further dilution, we recommend diluting these samples additionally 1:2 with REABUF.

Derivatization of standards (STD), controls (CTRL) and samples (SAMPLE) is carried out in single analysis in vials (e.g. 1.5 ml vials).

The reagents provided with this kit are sufficient for up to 48 derivatizations, which are transferred in duplicate determinations to the wells of the microtiter plate.

1. Bring all reagents and samples to **room temperature** (15-30 °C) and mix well.
2. Add **100 µl of diluted standards/ controls/ samples (STD/ CTRL/ SAMPLE)** in the corresponding vial.

3. Add **25 µl** of freshly prepared **derivatization reagent (DER)** into each vial (STD, CTRL, SAMPLE) and **mix thoroughly** by repeated inversion or several seconds on a vortex mixer. Incubate for **45 min at room temperature** (15-30°C) on a **shaker**.
4. Afterwards add **1000 µl of assay buffer (ASYBUF)** into each vial, mix well and incubate for **45 min at room temperature** (15-30°C) on a **shaker**.

2 x 25 µl of each treated sample (STD, CTRL, SAMPLE) are used in the ELISA as duplicates.

Test procedure

5. Mark the positions of standards/ controls/ samples (STD/ CTRL/ SAMPLE) in duplicate on a **protocol sheet**.
6. Take as many **microtiter strips** (PLATE) as needed from kit. Store unused strips covered at 2-8°C. Strips are stable until the expiry date stated on the label.
7. Wash each well **5 times** with **250 µl of wash buffer**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
8. For the analysis in duplicate, take **2 x 25 µl** of the **derivatized standards/ controls/ samples (STD/ CTRL/ SAMPLE)** out of the vials and add into the respective wells of the microtiter plate.
9. Add **100 µl of dissolved tyrosine antibody (AB)** into each well. Cover the plate tightly.
10. Incubate overnight (**15-20 hours**) at **2-8°C**.
11. Discard the content of each well and wash **5 times** with **250 µl of wash buffer**. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper.
12. Add **100 µl** of diluted **peroxidase labeled antibody (2. AB)** into each well.

13. Cover plate tightly and incubate for **1 hour at room temperature** (15-30°C) on a horizontal **shaker**.
14. Discard the content of each well and wash **5 times** with **250 µl of wash buffer**. After the final washing step the inverted microtiter plate should be firmly tapped on absorbent paper.
15. Add **100 µl** of **TMB substrate (SUB)** into each well.
16. Incubate for **8-12 min at room temperature** (15-30°C) in the dark*.
17. Add **100 µl of ELISA stop solution (STOP)** into each well, mix thoroughly.
18. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm (690 nm) as a reference.

* The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or Immundiagnostik AG.

7. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4 parameter algorithm".

1. 4 parameter algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e.g. 0.001).

2. Point-to-point calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before automatically evaluating the results. If this option is not available within the used program, the duplicate values should be evaluated manually.


EDTA plasma and serum

The concentrations can be determined directly from the standard curve in $\mu\text{mol/l}$. **No factor** is required.

The results from additionally diluted samples must be multiplied by this dilution factor.

In the following, an example of a calibration curve is given; do not use it for the calculation of your results.

Example of calibration curve

8. LIMITATIONS

Samples with an OD lower than the OD of the highest standard should be further diluted with reaction buffer (REABUF) and re-assayed. Please consider this dilution factor when calculating the results.

9. QUALITY CONTROL

Immundiagnostik recommends the use of external controls for internal quality control, if possible.

Control samples should be analyzed with each run. Results generated from the analysis of control samples should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control samples are outside of the acceptable limits.

Reference Range

Based on internal studies with plasma samples of apparently healthy persons (n=146) a mean value of 58 µmol/l was calculated. The standard deviation was 14.4 µmol/l.

Plasma mean value ± 2 × standard deviation: **58 ± 28.8 µmol/l**

Normal range: **29 – 87 µmol/l**

We recommend each laboratory to establish its own reference range.

10. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-assay (n=19)

Sample	Tyrosine [µmol/l]	CV [%]
1	31.4	7,6
2	125.3	7.4

Inter-assay (n=8)

Sample	Tyrosine [µmol/l]	CV [%]
1	69.7	5.8
2	104.6	4.6

Spiking recovery

Two plasma samples were spiked with different tyrosine concentrations and measured in this assay. The mean recovery rate for all concentrations was 93.2 % (n=5).

Sample [µmol/l]	Spike [µmol/l]	Tyrosine expected [µmol/l]	Tyrosine measured [µmol/l]	Recovery [%]
76.8	25	101.8	94.5	92.8
76.8	50	126.8	123.3	97.2
108.8	25	133.8	116.1	86.8
108.8	50	158.8	152.4	96.0

Dilution recovery

Two spiked plasma sample was diluted and measured in this assay. The mean recovery was 97.1 % (n=5).

Sample [µmol/l]	Dilution	Expected [µmol/l]	Measured [µmol/l]	Recovery [%]
76.8	1:2	38.4	34.0	88.4
76.8	1:4	19.2	20.0	103.9
108.8	1:2	54.4	50.6	93.0
108.8	1:4	27.2	28.0	102.9

Analytical sensitivity

The zero-standard was measured 82 times. The detection limit was set as $B_0 - 2 \text{ SD}$ and estimated to be 3.6 µmol/l.

Sample	Tyrosine mean value [OD]	2 x standard deviation (SD)	Detection limit [µmol/l]
Zero-standard	1.24	0.15	3.6

Specificity

Specificity was tested by measuring the cross-reactivity against compounds with structural similarity to tyrosine. The specificity is calculated in percent in relation to the tyrosine binding activity.

L-phenylalanine	< 2 %
L-tryptophan	< 0.5 %

11. PRECAUTIONS

- All reagents in the kit package are for research use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or ProClin as bactericides. Sodium azide and ProClin are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of sulfuric acid, which is a strong acid. Although diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped up immediately with copious quantities of water. Do not breathe vapour and avoid inhalation.

12. TECHNICAL HINTS


- Do not interchange different lot numbers of any kit component within the same assay. Furthermore, we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch, as wells from already opened microtiter plates are exposed to different conditions than sealed ones.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results proper adhesion of plate sealers during incubation steps is necessary.

- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature, and pipetting volumes of the different components are defined by the producer. Any variation of the test procedure which is not coordinated with the producer may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from incorrect use.
- Warranty claims and complaints regarding deficiencies must be logged within 14 days after receipt of the product. The product should be sent to Immundiagnostik AG along with a written complaint.

Used symbols:

	Temperature limitation		Catalogue Number
	For research use only		To be used with
	Manufacturer		Contains sufficient for <n> tests
	Lot number		Use by
	Attention		