ClinCal® - Calibrator
Plasma Calibrator, lyophilised, for Mycophenolic Acid
Plasma-Kalibrator, lyophilisiert, für Mycophenolsäure

for HPLC Assay in Plasma

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REF	29013
LOT	1186
	2020-05

Analyte / Analyt	Mean value / Sollwert
Mycophenolic Acid / Mycophenolsäure	4.29 mg/l (13.4 µmol/l)

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Intended use:

ClinCal® Plasma Calibrators are used for calibration of the HPLC-system. After reconstitution this lyophilised calibrator has to be prepared like a patient sample.

Reconstitution:

Add exactly 3.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma calibrator is stable for 48 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analyte is:

at least 8 hours when stored at 15 - 30 °C
 at least 7 days when stored at 2 - 8 °C
 at least 30 days when stored below -18 °C
 (avoid repeated freezing and thawing)

The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.87 %.

Pack size:

ClinCal® Plasma Calibrator for Mycophenolic Acid: 1 x 3 ml, **order no.: 29013**

Caution:

The human plasma which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the plasma calibrator should be considered as potentially infectious and treated with appropriate care.

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