## ClinCal® - Calibrator Whole Blood Calibrator lyophilised / Vollblut-Kalibrator lyophilisiert

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### Intended use:

ClinCal® Whole Blood Calibrators are used for calibration of the analysis-system. This Level 7 shall be used in combination with the Calibrator Set, order no. 9933, or the Calibrator Set, order no. 9033, for expansion of the calibration range. After reconstitution the calibrator has to be prepared like patient samples in one series of analyses.

### Reconstitution:

Add exactly 2.0 ml of HPLC-water to the vial. To dissolve the material completely, we recommend a so-called roller mixer. A treatment of 30 to 60 minutes on such an instrument is sufficient.

### Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 30 days when stored below -18 °C (avoid repeated freezing and thawing)

### Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over an extended range of concentration. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.56 %

### Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK) with statistical methods. The mean values for Tacrolimus and Sirolimus are traceable to ERM®-DA110a and ERM®-AC021a, respectively. Those for Cyclosporine A and Everolimus to CRM from Cerilliant Corp.

### Pack size:

ClinCal® Whole Blood Calibrator for Immunosuppressants Level 7  
2 x 2 ml, **order no.: 9028**

### Precautions:

The human whole blood which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the calibrator should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH