ClinChek® - Control

Urine Control lyophilised / Kontrollurin lyophilisiert

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Intended use:

ClinChek® urine controls are used for internal quality assurance in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. These lyophilised controls are based on human urine and are available in three different levels of concentration. After reconstitution the controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 3.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised urine control is stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 12 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.98 %.

The creatinine concentration of the different control levels is 0.23 g/l.

Mean values:

The mean values and confidence intervals have been established at RECIPE, according to the guideline of the German Medical Association for quality assurance (Rili-BÄK), with statistical methods.

Pack size:

ClinChek® Urine Control
Level I
5 x 3 ml, **order no.: MS3080**
Level II
5 x 3 ml, **order no.: MS3081**
Level III
5 x 3 ml, **order no.: MS3082**
Level I, II, III
3 x 2 x 3 ml, **order no.: MS3083**

Precautions:

These controls were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH