ClinCal[®] - Calibrator

Whole Blood Calibrator Set lyophilised / Vollblut-Kalibrator-Set lyophilisiert
FOR IMMUNOSUPPRESSANTS / FÜR IMMUNSUPPRESSIVA

REF	9933
LOT	2430
	2023-10

Intended use:

ClinCal[®] Whole Blood Calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human blood and are available as a calibrator set with 7 levels, including a blank. After reconstitution the calibrators have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 2.0 ml of HPLC-water to the vial. To dissolve the material completely, we recommend a so-called roller mixer. A treatment of 30 to 60 minutes on such an instrument is sufficient.

Storage and stability:

This product will be stable until the expiration date when stored unopened at 2 - 8 °C.

After reconstitution the stability of the analytes is:

- 1 day when stored at 15 - 30 °C
- 7 days when stored at 2 - 8 °C
- 28 days when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a wide range of concentration. According to quality assurance all ClinCal[®] Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.61 %

Mean values:

The mean values have been established at RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. The mean values for Tacrolimus and Sirolimus are traceable to ERM[®]-DA110a and ERM[®]-DA111a, respectively. Those for Cyclosporine A and Everolimus to Cerilliant CRM C-093 and E-068, respectively.

Pack size:

ClinCal[®] Whole Blood Calibrator Set for Immunosuppressants
Level 0 - 6
7 x 1 x 2 ml, **order no.: 9933**

Precautions:

The human whole blood which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV1/2-, HBV- and HCV-antibodies, Hepatitis B-surface antigen, HIV1- and HCV-RNA, HBV-DNA (NAT). Nevertheless, the calibrators should be considered as potentially infectious and treated with appropriate care.

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Level 0

This blank-level was manufactured from the same pool of human blood like the others, but it was not spiked with any Immunosuppressant.

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In Germany Cyclosporine A and Tacrolimus are subject to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK). For assessment of the internal quality control, please also refer to table B1a of the guideline.

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