ClinCal® - Calibrator

Plasma Calibrator lyophilised / Plasma-Kalibrator lyophilisiert
FOR TRACE ELEMENTS / FÜR SPURENELEMENTE

REF	9985
LOT	1177
	2022-04

Intended use:

ClinCal® Plasma Calibrators are used for calibration of your analyzing system in the occupational and environmental toxicological laboratory as well as in the clinical-chemical one. This lyophilised calibrator is based on human plasma and is available in a range of concentration that allows a comfortable one-point-calibration. After reconstitution the calibrator has to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 3.0 ml of analytical grade or double distilled water to the vial and close it well. Then dissolve the lyophilised material completely. Prior to use mix the calibrator well again.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma calibrator is stable for 48 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30°C
- at least 4 days when stored at 2 - 8 °C
- at least 1 month when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The concentrations of the analytes are chosen in ranges where valid results can be obtained. According to quality assurance all ClinCal® Calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume (CV) is < 1 %. The average residual moisture of this lot is 0.98 %.

Mean values:

The mean values have been established in independent reference laboratories with supervision of RECIPE, according to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK), with statistical methods. For the determination of the trace elements, atomic absorption spectrometry (AAS) was used as well as ICP-MS. Even if you use different techniques of analysis, you should find your control results within the ranges indicated on the data sheet, after having calibrated your analyzing system with this ClinCal® plasma calibrator.

Pack size:

ClinCal® Plasma Calibrator
5 x 3 ml, **order no.: 9985**

Precautions:

The human plasma which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV1/2- and HCV-antibodies, Hepatitis B-surface antigen. Nevertheless, the plasma calibrator should be considered as potentially infectious and treated with appropriate care.

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For the following elements we could not establish any mean values, due to the low number and the high imprecision of the available results: **Au, Be, Fe, Fluoride, Li, Mo, Ti.**

For more information please contact us : customerservice@recipe.de /

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In Germany Lithium and Magnesium are subject to the Guideline of the German Medical Association on Quality Assurance (Rili-BAEK). For assessment of the internal quality control, please also refer to table B 1a of this guideline. /

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